





Blood Groups 🔍



The ABO blood group system classifies blood into four types (A, B, AB, O) based on specific antigens on red blood cells, while the Rh factor (positive or negative) further classifies blood depending on the presence of the Rh antigen.

- * Blood group A has A antigens on the red blood cells with anti-B antibodies in the plasma
- * Blood group B has B antigens on the red blood cells with anti-A antibodies in the plasma
- * Blood group O has no antigens, but both anti-A and anti-B antibodies in the plasma
- * Blood group AB has both A and B antigens, but no antibodies



CUSTOMER FEATURES 🧈





* Allows new users to create an account by providing basic information for secure access to the application.



* Enables users to securely log in with their credentials to access blood group identification features.



* Serves as the main dashboard, offering a simple interface to navigate blood identification tools and user options.



- * Displays user information and blood group details.
- * Allows users to upload images of ABO blood cells to check the blood group.

MILESTONES The system allows secure user The extracted contour count is registration, login, and image displayed on the profile page, uploads for blood group integrated into the blood group identification. It also stores and identification process. This provides manages user credentials in a users with accessible and relevant structured database. information for their results. 01 03 02 04 The system uses OpenCV to The system preprocesses blood preprocess blood slide images with slide images with OpenCV and grayscale conversion, Gaussian applies morphological operations blur, and thresholding to enhance to identify the ABO blood type features. It then detects contours and Rh factor. The processed


to analyze key elements for blood


group identification.

image and results are displayed

on the profile page for the user.




















Advantages Provides faster, more accurate blood typing while eliminating manual errors. Delivers quick results suitable for both emergency and routine cases. 4 Utilizes OpenCV for precise

Designed for accessibility by both medical professionals and regular users without requiring specialized knowledge. Utilizes OpenCV for precise image processing and is cost-effective, scalable, and web-based for secure remote

The Blood Group Identification System provides an efficient, automated solution for identifying blood groups through blood cell image uploads. Leveraging OpenCV with techniques like grayscale conversion, thresholding, and morphological operations, it ensures accurate ABO and Rh factor detection. The user-friendly web interface simplifies blood group identification, making it accessible to both medical professionals and individuals while

reducing manual errors.



